
**DE** Gebrauchsanweisung **EN** Instructions for Use **FR** Mode d'emploi  
**NL** Gebruiksaanwijzing **IT** Istruzioni d'uso **TR** Kullanma Kılavuzu  
**ES** Manual de instrucciones **RU** Инструкция по пользованию **EL** Οδηγίες χρήσης  
**SK** Návod na použitie **DA** Brugsanvisning **SV** Bruksanvisning **ZH** 使用说明书










# NP15

Nasal Pillow Mask

# NP15

Nasal Pillow Mask

<b>Deutsch</b>	<b>3</b>
<b>English</b>	<b>11</b>
<b>Français</b>	<b>19</b>
<b>Nederlands</b>	<b>27</b>
<b>Italiano</b>	<b>35</b>
<b>Svenska</b>	<b>44</b>
<b>Русский</b>	<b>51</b>
<b>Türkçe</b>	<b>60</b>
<b>Ελληνικά</b>	<b>68</b>
<b>中文</b>	<b>76</b>
<b>Español</b>	<b>82</b>
<b>Dansk</b>	<b>91</b>
<b>Slovensky</b>	<b>99</b>

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---

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## Contents

---

<b>1 Overview</b> .....	<b>11</b>
2.1 Safety instructions .....	11
2.2 Contraindications .....	12
2.3 Side effects .....	12
<b>3 Description of device</b> .....	<b>12</b>
3.1 Intended use .....	12
3.2 Duty to provide instruction .....	12
3.3 Description of function .....	12
<b>4 Operation</b> .....	<b>13</b>
4.1 Put on mask .....	13
4.2 Remove mask .....	14
4.3 Dismantle/assemble mask .....	14
<b>5 Hygiene treatment</b> .....	<b>14</b>
5.1 Procedure .....	14
5.2 Permitted processes .....	15
5.3 Description of processes .....	15
5.4 Following treatment .....	15
5.5 Change of patient .....	15
<b>6 Function check</b> .....	<b>16</b>
6.1 Intervals .....	16
6.2 Procedure .....	16
<b>7 Servicing</b> .....	<b>16</b>
7.1 Storage .....	16
7.2 Disposal .....	16
<b>8 Troubleshooting</b> .....	<b>16</b>
<b>9 Technical data</b> .....	<b>16</b>
9.1 Specifications .....	16
9.2 Materials .....	17
9.3 Characteristic pressure/flow curve ..	17
<b>10 Warranty</b> .....	<b>17</b>
<b>11 Declaration of conformity</b> .....	<b>18</b>

## 1 Overview

---

For individual parts, see enclosed leaflet (Fig. A):

- 1 Headgear
- 2 Crossband
- 3 Tube anchoring clip
- 4 Forehead support
- 5 Mask body
- 6 Elbow
- 7 Tube adapter
- 8 Pressure measurement port with plug
- 9 Tube anchoring strap
- 10 Mask tube
- 11 Rotating sleeve
- 12 Headgear clip
- 13 Nasal pillow
- 14 Forehead cushion
- 15 Forehead strap
- 16 Cheek strap

## 2 Safety instructions

---

### 2.1 Safety instructions

Read these instructions for use through carefully. They are a constituent part of the mask and must be available at all times. In accordance with the requirements of directive 93/42/EEC, observe the points below.

#### 2.1.1 Hazards

##### **Risk of injury from additionally-supplied oxygen!**

Oxygen can become deposited in clothing, bed linen and hair. In conjunction with smoking, naked flame and electrical equipment, it can cause fires and explosions.

- Do not smoke.

- Do not use naked flame.
- Do not use electrically-conductive tubes.

### **Risk of injury from oxygen therapy!**

Oxygen therapies can lead to side effects if incorrectly used/metered.

- Only perform oxygen therapy following prescription by a physician.
- Follow the instructions for use of the oxygen system and the therapy device.

### **Risk of injury from re-inhaling CO<sub>2</sub>!**

If the mask is used incorrectly, CO<sub>2</sub> may be re-inhaled.

- Only put mask on for a prolonged period if the therapy device is switched on.

### **Risk of injury from mask slipping!**

If the mask slips or falls off, therapy is ineffective.

- Subject patients with restricted ability to breathe spontaneously to continuous monitoring.
- Activate low pressure/leakage alarms on ventilation device.
- Subject patients with obstructive and restrictive lung diseases to individual monitoring.

#### **2.1.2 General**

- If third-party items are used, functional failures and restricted fitness for use may result. Biocompatibility requirements may also not be met. Please note that in these cases, any claim under warranty and liability will be void if neither the accessories nor genuine replacement parts recommended in the instructions for use are used.
- To prevent an infection, bacterial contamination or functional impairments, follow the instructions in the section entitled "Hygiene treatment" on Page 14.

## **2.2 Contraindications**

The following contraindications are possible: serous rhinitis, nosebleeds, erosions and ulcerations around the nostrils, skin allergies.

## **2.3 Side effects**

Blocked nose, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated mucous membrane, skin rashes, irritating noises when breathing.

## **3 Description of device**

---

### **3.1 Intended use**

Nasal pillow mask NP15 forms a link between the patient and the therapy device. It is used to treat sleep apnea and patients with ventilatory insufficiency. The mask is not intended for life-support ventilation.

### **3.2 Duty to provide instruction**

Provide the patient with instruction in how to use the mask before therapy starts.

### **3.3 Description of function**

#### **3.3.1 Therapy pressure range**

The therapy pressure range is 4 hPa - 25 hPa.

#### **3.3.2 Compatible devices**

- You can use the mask with all therapy devices which are not for life-support treatment of patients with sleep apnea or ventilatory insufficiency.
- With many combinations of device, the actual pressure in your mask will not correspond to the prescribed therapy pressure, even if the therapy device is displaying the correct pressure. Have the device combination adjusted by a physician or specialist dealer so that the actual pressure in the mask corresponds to therapy pressure.

- If you are using the mask with a Weinmann sleep therapy device with an Autostart function, the therapy device may not always switch off automatically when the mask is taken off. This applies to size S nasal pillows in particular.
- If your ventilation device has leakage or low pressure alarms, set these so that an alarm is triggered as soon as you remove the mask or it slips out of position.

### 3.3.3 Exhalation system

The mask has an integrated exhalation system. Rotating sleeve **11** and elbow **6** are shaped so that there is a gap between these parts (Fig. **G**), through which exhaled air can escape. The specific geometry of the parts ensures that the rotating sleeve cannot seal this gap when it is pushed into place.

### 3.3.4 Pressure measurement port

The mask has a pressure measurement port **8** on tube adapter **7** with which you can either measure therapy pressure or supply oxygen. If you are not using the pressure measurement port, seal it with the plug to reach therapy pressure.

### 3.3.5 Respiratory air humidifier

The mask can be used with cold or warm air humidifiers, though the degree of humidification should never be so high that water condenses in the mask tube. Follow the instructions for use of the respiratory air humidifier.

## 4 Operation

All the figures are on the enclosed leaflet.

### 4.1 Put on mask

1. Select a suitable nasal pillow.
2. Push nasal pillow onto mask body (Fig. **B**).
3. Push nasal pillow against nostrils.
4. Pull headgear over head (Fig. **K**).
5. Using the Velcro tab on crossband **2** on the head, adjust the width of the headgear so that the crossband is located centrally on the head (Figs. **H-J**).
6. Adjust forehead strap **15** and cheek strap **16** using the Velcro tabs on both sides so that the mask seals on the face, but is not too tight.
7. Keep the tabs of forehead support **4** depressed.
8. Engage the forehead support in one of the seven positions (numbers 1-7, Fig. **C**) so that the mask seals on the face but is not too tight.
9. Connect mask tube **10** to the ventilation tube of the therapy device via tube adapter **7**.

#### Tip

- The ball joint allows you to guide the mask tube freely (Fig. **I**). You can attach it to the headgear on the side (Fig. **J**) or top (Fig. **H**) using the tube anchoring strap.
- If you wish to fix the mask tube to the head for sleeping, first connect the mask tube to the ventilation tube of the therapy device. You can then either attach the mask tube to the forehead support using tube anchoring clip **3** and/or to the top of the headgear using tube anchoring strap **9**.
- Ensure that there is a secure, airtight connection to the ventilation tube. The connection may not come loose as a result of therapy.

## 4.2 Remove mask

1. Remove mask without changing settings.
2. Disconnect mask tube from the ventilation tube of the therapy device.

## 4.3 Dismantle/assemble mask

1. Take plug off pressure measurement port (Fig. **F**).
2. Release mask tube and rotating sleeve **11** from elbow.
3. Release nasal pillow from mask body.
4. Release forehead cushion from forehead support by pulling sideways (Fig. **D**).
5. Release headgear from headgear clips and forehead support.
6. Subject mask parts and headgear to a hygiene treatment.
7. Following the hygiene treatment, reassemble the mask in reverse sequence.


### Tip

- The Velcro tabs on forehead strap **15** are pale blue, the Velcro tabs on cheek strap **16** are dark blue (Fig. **E**).
- Always thread the forehead strap for the headgear through slot a (Fig. **Ea**) of the forehead support and the tube anchoring clip through slot b (Fig. **Eb**) of the forehead support.

## 5 Hygiene treatment

---

### 5.1 Procedure

#### Risk of injury if cleaning is inadequate!

Residues may block mask and mask tube. They may impair the integrated exhalation system and jeopardize the success of therapy.

- For patients with a compromised immune system or particular background of illness, disinfect or sterilize mask parts daily following consultation with the physician.
- Carefully clean the contact surfaces between the elbow and the rotating sleeve.

#### 5.1.1 Daily

1. Dismantle mask as described.
2. Wash mask in hot water using a mild detergent (e.g. washing-up liquid).
3. Carefully rinse off all parts with clean water.  
The parts must be free of residues.
4. Allow parts to air-dry.

#### 5.1.2 Weekly

1. Dismantle mask as described.
2. Subject mask and headgear to a hygiene treatment (see table).

### Tip

- Wear suitable gloves (e.g. disposable gloves) during the hygiene treatment.
- Do not dry parts in direct sunlight.
- Wash the headgear before using for the first time, as the dye may run.
- Do not iron the headgear, otherwise the Velcro tabs will stop gripping.
- Do not dry the headgear in a dryer.

## 5.2 Permitted processes

Process	Dishwasher	Hand wash	Disinfect <sup>1</sup>	Sterilize <sup>2</sup>
Nasal pillow	•	•	•	•
Forehead cushion	•	•	•	•
Plug	•	•	•	•
Mask body	•	•	•	
Elbow	•	•	•	
Headgear clip	•	•	•	
Headgear		•		
Mask tube	•	•	•	
Tube anchoring clip	•	•	•	
Tube anchoring strap		•		

<sup>1</sup> Compatibility has been proven for the disinfectants GIGASEPT® FF, CIDEX® OPA and ANIOSYME® DD1. Follow the instructions for use for the disinfectants.

<sup>2</sup> Alternative to disinfecting.

## 5.3 Description of processes

Dishwasher	Cycle at 65 °C. Then carefully rinse off all parts in clean water.
Hand wash	Hand wash in warm water. Use a mild detergent. Then carefully rinse off all parts with clean water. Allow to air-dry.
Disinfect <sup>1</sup>	Clean mechanically in dilute solution, using a brush, for example. Carefully rinse off all parts with clean water. Allow to air-dry.
Sterilize (alternative to disinfecting)	Steam-sterilize in devices to EN 285. Temperature 134 °C, minimum dwell time 3 minutes.

Notes on the hygiene treatment and on the hygiene treatment for a change of patient can be found in a brochure from our information portal at [www.loewensteinmedical.de](http://www.loewensteinmedical.de). If you do not have Internet access, contact Löwenstein Medical.

## 5.4 Following treatment

1. Thoroughly dry all parts following the hygiene treatment.
2. Repeat the hygiene treatment if required.
3. Perform a function check.
4. Reassemble mask.

## 5.5 Change of patient

1. Dismantle mask.
2. Wash mask in hot water using a mild detergent (e.g. washing-up liquid) so that no residues remain before disinfecting or sterilizing.
3. Subject mask to a hygiene treatment (see table).
4. Use new headgear.

## 6 Function check

### 6.1 Intervals

Perform a function check after every hygiene treatment.

### 6.2 Procedure

If the forehead cushion becomes discolored, this does not affect its function.

1. Check all parts for soiling and damage. Replace parts if necessary.
2. In the event of leaks due to damage, replace the nasal pillow.

## 7 Servicing

The mask requires no servicing.

### 7.1 Storage

The mask may only be stored under the ambient conditions described (see "Technical data" on Page 16)

### 7.2 Disposal

Parts which are no longer usable can be disposed of with domestic waste.

## 8 Troubleshooting

Fault	Cause	Remedy
Pain in the face or on the nose from pressure	Mask too tight	Adjust headgear and modify position of forehead support

Fault	Cause	Remedy
Therapy pressure not reached	Mask not correctly adjusted	Adjust headgear and modify position of forehead support
	Air escaping between nasal pillow and nose	Position nasal pillow correctly or use nasal pillow of a different size
	Nasal pillow damaged	Replace nasal pillow
	Tube system leaking	Check connector and firm location of tube sleeves
	Air escaping at pressure measurement port	Seal pressure measurement port with plug

## 9 Technical data

### 9.1 Specifications

NP15	Value
Product class to directive 93/42/EEC:	II a
Dimensions (W x H x D):	approx. 140 mm x 100 mm x 40 mm
Weight:	approx. 90 g <sup>1</sup>
Therapy pressure range:	4 hPa - 25 hPa
Tube connection: tapered connector to DIN EN ISO 5356-1	Ø 22 mm (male)
Pressure measurement port:	Ø 4 mm
Temperature range Operation: Storage:	+ 5 °C to + 40 °C - 20 °C to + 70 °C

NP15	Value
Flow resistance <sup>3</sup> at 50 l/min at 100 l/min	≤ 0.9 hPa ≤ 3.5 hPa
Quoted two-number noise emission value to ISO 4871: – sound pressure level – sound power level (uncertainty factor: 3 dB(A))	15 dB(A) 23 dB(A)
Service life	up to 12 months <sup>2</sup>

<sup>1</sup> Depending on size

<sup>2</sup> Materials used to manufacture masks will age if exposed to aggressive detergents, for example. In individual cases, therefore, it may be necessary to replace your mask sooner (see "6. Function check" on Page 16).

<sup>3</sup> Flow resistance depends on the size of the nasal pillow and may deviate.


## 9.2 Materials

Article	Material
Nasal pillow, pressure measurement port plug, forehead cushion	Silicone
Mask body, rotating sleeve, tube adapter (ball socket), forehead support, headgear clips, Velcro tabs	Polyamide
Elbow, tube adapter (ball)	Polycarbonate
Mask tube	Ethylene vinyl acetate copolymer
Headgear, tube anchoring strap	Lycra, polyester, polyurethane, UBL loop

All parts of the mask are free from latex, PVC and DEHP.

## 9.3 Characteristic pressure/flow curve

The characteristic pressure/flow curve shows outlet flow rate as a function of therapy pressure.

**CE 0197** The right to make design modifications is reserved.

## 10 Warranty

- Löwenstein Medical warrants that the product will remain free from defects for a period of six months from the date of purchase if used in accordance with purpose.
- A condition for claims under warranty is presentation of a proof of purchase showing salesperson and date of purchase.
- We do not accept claims under warranty if:
  - the instructions for use are not followed
  - there are operating errors
  - the product is not used or treated properly
  - an unauthorized third party interferes with the device for repair purposes
  - there is force majeure
  - transport damage is sustained due to incorrect packing of returned goods
  - normal wear and tear is sustained during operation
  - genuine replacement parts are not used.
- Löwenstein Medical accepts no liability for consequential damages due to

defects unless these are due to deliberate or gross negligence or in the case of negligent injury to life or limb.

- Löwenstein Medical reserves the right to decide whether to rectify the defect, supply a defect-free item or reduce the purchase price by a reasonable amount.
- If we reject a claim under warranty, we will not bear any of the expense of transport between customer and manufacturer.
- Statutory claims under warranty are not affected by these provisions.

## **11 Declaration of conformity**

---

Löwenstein Medical GmbH + Co. KG hereby declares that the product complies with the relevant provisions of directive 93/42/EEC pertaining to medical devices.

The complete text of the declaration of conformity can be found at:  
[www.loewensteinmedical.de](http://www.loewensteinmedical.de)

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