

Component	Material
Mask Base	Polycarbonate (PC)
Non-Rebreathing Valve	Polycarbonate (PC)/Silicone (LSR)
Elbow	Polycarbonate (PC)
Glider Strap	Nylon (PA)
Glider Clip	Acetal (POM)
Forehead Cushion	Silicone (LSR)
Silicone Seal (S, M, L, XL)	Silicone (LSR)
Foam Cushion (S, M, L, XL)	Polyurethane
Headgear	Nylon/Polyurethane
Pressure Ports Connector (Euro Only)	Polycarbonate (PC)/Santoprene (TPE)