

Quattro™ Air

Full Face Maske

Komponente	Material	
Maskenrahmen	Hart	Nylon
	Weich	Silikon-Elastomer
Kniestück	Hart (transparent)	Nylon
	Hart (blickdicht)	Polypropylen
	Weich	Thermoplastisches-Elastomer
	Ventil	Silikon-Elastomer
Maskenkissen	Hart	Nylon
	Weich	Silikon-Elastomer
Kopfband	Textil	Nylon/Lycra
	Flexible Bänderung	Polyurethanschaum